

Ортодонтические инструменты HSL, HSJ, HSG, HSC, HWC, HWL, HWM, HSK, HSB, HSE

Описание

По вопросам продаж и поддержки обращайтесь:

Алматы (727)345-47-04
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(727)345-47-04

Беларусь +(375)257-127-884

Узбекистан +998(71)205-18-59


Киргизия +996(312)96-26-47

эл.почта: hrg@nt-rt.ru || сайт: <https://hammacher.nt-rt.ru/>

HSL 960-02, Stand for pliers

Category: Stands for pliers

Order number: HSL_960-02




Hünebeck

HSL 958-03, Stand for pliers

Category: Stands for pliers

Order number: HSL_958-03




Hünebeck

HSL 958-02, Stand for pliers

Category: Stands for pliers

Order number: HSL_958-02




Hünebeck

HSL 958-01, Stand for pliers

Category: Stands for pliers

Order number: HSL_958-01



Hünebeck

HSL 497-14, Side Nippers

Category: Pliers

Order number: HSL_497-14




disinfectable 140 mm / 5-1/2 inch chromium plated For hard wire until max. 1,2 mm / .047" diameter

HSL 483-13, Flat nose pliers

Category: Pliers

Order number: HSL_483-13



sterilizable with guide notch, serrated 130 mm / 5-1/8 inch For hard wire until max. 0,7 mm / .028" max.

HSL 481-12, Flat nose pliers

Category: Pliers

Order number: HSL_481-12



sterilizable narrow, serrated 120 mm / 4-3/4 inch For hard wire until 0,5 mm / .020" diameter max.

HSL 478-13, Flat nose pliers

Category: Pliers

Order number: HSL_478-13



sterilizable stainless, serrated 135 mm / 5-3/8 inch For Wire until 0,7 mm / .028" max.

HSL 446-12, Tweed arch pliers

Category: Pliers

Order number: HSL_446-12



Angle/Tweed, sterilizable 120 mm / 4-3/4 inch 2 steps Bends round, square and rectangular wires up to 0,56 mm x 0,7 mm / .020 x .028"

HSL 436-15, Universal pin bending pliers

Category: Pliers

Order number: HSL_436-15



sterilizable 155 mm / 6-1/8 inch with very delicate working ends Bends: hard wire up to 0,9 mm / .035" diameter max. Cuts: hard wire up to 0,7 mm / .028" diameter max.

HSL 431-15, Universal pliers

Category: Pliers

Order number: HSL_431-15



Waldsachs, sterilizable 150 mm / 5-7/8 inch Bends: hard wire up to 1,0 mm / .039" diameter max. hard profiled wires up to 1,3 x 0,65 mm / .051" x .026" max. Cuts: hard wire up to 0,7 mm / .028" diameter max.

HSL 4292-13, Arch contouring pliers

Category: Pliers

Order number: HSL_4292-13



De La Rosa arch contouring pliers, non-grooved, sterilizable 120 mm / 4-3/4 inch Non-grooved pliers allow countouring the archwire at any position on the instrument. Bends round, square and rectangular wires up to 0,5 mm / .020".

HSL 4291-13, Arch contouring pliers

Category: Pliers

Order number: HSL_4291-13



De La Rosa arch contouring pliers, grooved, sterilizable 120 mm / 4-3/4 inch Precision machined grooves measure 0,41 mm 0,46 mm and 0,56 mm .016" .018" and .022".

HSL 4247-13B, Torquing pliers

Category: Pliers

Order number: HSL_4247-13B



125 mm (4 3/4 inch), No. d261B, with TC-inserts

HSL 4247-13A, Torquing pliers

Category: Pliers

Order number: HSL_4247-13A



125 mm (4 3/4 inch), No. d261A, with TC-inserts

HSL 4247-13, Torquing pliers

Category: Pliers

Order number: HSL_4247-13




Set from HSL 4247-14A and HSL 4247-14B, 125 mm / 4-7/8, sterilizable

HSL 4246-13, Torque-Set

Category: Pliers

Order number: HSL_4246-13



Single tooth torque set for anterior mandibles for 3,5 mm bracket width.125 mm / 4-7/8 inch, sterilizableThe othodontist can set the required torque directly during the removal of the arch without twisting the remaining wire.A displacement resp. debonding

HSL 4243-13, Step pliers

Category: Pliers

Order number: HSL_4243-13



sterilizable With TC-inserts 125 mm / 4-7/8 inch Automatic bending of 0,75 mm steps, For round, square and rectangular wire until max. 0,53 x 0,64 mm / .021 x .025"

HSL 4240-13, Stop pliers V-bend

Category: Pliers

Order number: HSL_4240-13



sterilizable With TC-inserts 125 mm / 4-7/8 inch For placing a 1 mm stop For round, square and rectangular wire until max. 0,56 mm / .022"

HSL 4192-12B Wire and clasp bending pliers

Category: Pliers

Order number: HSL_4192-12_B



Hammacher / Aderer, sterilizable 120 mm / 4-3/4 inch 90° angled For hard wire until 0,9 mm / .035" diameter max.

HSL 4192-12 A Wire and clasp bending pliers

Category: Pliers

Order number: HSL_4192-12A



Hammacher / Aderer, sterilizable 120 mm / 4-3/4 inch 90° angled For hard wire until 0,9 mm / .035° diameter max.

HSL 4191-12B Aderer

Category: Pliers

Order number: HSL_4191-12B



Wire and clasp bending pliers, 90° angled For hard wire until 0,7 mm /.028" diameter

HSL 404-15, Sling, clasps and wire bending pliers

Category: Pliers

Order number: HSL_404-13



Adams, sterilizable 130 mm / 5-1/8 inch For hard wire up to 0,7 mm / .028" diameter max.

HSL 403-12, Sling, clasps and wire bending pliers

Category: Pliers

Order number: HSL_403-12



Angle modif., sterilizable 125 mm / 4-7/8 inch For hard wire up to 0,7 mm /.028" diameter

HSL 4015-13, Wire bending pliers

Category: Pliers

Order number: HSL_4015-13



Angle/Bird Beak, serrated, sterilizable With TC-inserts 125 mm / 4-7/8 inch For bending of round, square and rectangular wire up to max. 0,76mm / .030".

HSL 4011-13, Wire bending pliers

Category: Pliers

Order number: HSL_4011-13



Angle/Bird Beak, no. 139, sterilizable 120 mm / 4-3/4 inch Bends round, rectangular and square wires up to 0,76 mm /.030" Excellent for all types of wire bending in the laboratory. Stubby beaks - round on square.

HSL 3881-14, Light wire bending pliers

Category: Pliers

Order number: HSL_3881-14



Jarabak, sterilizable 130 mm / 5-1/8 inch Bends round, square and rectangular wires up to 0,5 mm / .020". Fine serrations on square beak and three precision grooves allow accurate bending & closing of loops with a secure grip. Especially good when used on

HSL 3850-13 Wire arches bending pliers "Nice End Mini"

Category: Pliers

Order number: HSL_3850-13



Nice End Mini, sterilizable 140 mm / 5-1/2 inch Suitable for wire cross-sections until max. 0,4mm / .016"

HSL 385-13, Wire arches bending pliers "Nice End"

Category: Pliers

Order number: HSL_385-13



Nice Ende, sterilizable 140 mm / 5-1/2 inch Suitable for wire cross-section until max. 0,46 mm / .018"

HSL 380-13, Optical pliers

Category: Pliers

Order number: HSL_380-13



Schwarz, sterilizable 130 mm / 5-1/8 inch For hard wire up to 0,7 mm / .028" diameter max.

HSL 369-13, Wire bending pliers

Category: Pliers

Order number: HSL_369-13



Adams mod., sterilizable 130 mm / 5-1/8 inch For hard wire until max. 0,9 mm / .035" diameter With one groove in front of the tip. Especially for triangular clasps, Adams clasps and springs of the protrusion.

HSL 364-13, Wire and bar bending pliers

Category: Pliers

Order number: HSL_364-13



Jarinda, sterilizable 125 mm / 4-7/8 inch For all commercially available orthodontic wires until max. 0,8 mm / .031" (hard and spring hard) Universally applicable, allows easy and comfortable working in the laboratory

HSL 344-00, Replacement tip made of nylon

Category: Pliers

Order number: HSL_344-00



For HSL 221-12, sterilizable 1 Piece

HSL 336-99, Replacement tip

Category: Pliers

Order number: HSL_336-99



For HSL 336-13, sterilizable 1 Piece

HSL 336-13, Band removing pliers

Category: Pliers

Order number: HSL_336-13



135 mm / 5-3/8 inch sterilizable with replacable tips

HSL 334-17, Bracket removing pliers

Category: Pliers

Order number: HSL_334-17



These pliers are especially appropriate for the lingual technique. Bracket removing pliers acc. to Dr. Wilhelmy, 170 mm / 6-3/4 inch Sterilizable

HSL 2842-13, Hard wire cutter

Category: Pliers

Order number: HSL_2842-13



With TC-inserts, 120 mm (4 3/4 inch). For hard wire from min. 0,6 mm / .024" until max. 0,8 mm / .031" diameter. With a 15 °angle . For an easier intra-oral access to the posterior area.

HSL 2724-13, Wire seizing and ligature pliers

Category: Pliers

Order number: HSL_2724-13



Weingart, sterilizable 135 mm / 5-3/8 inch For holding of wires round, square and rectangular until max. 0,56 mm x 0,70 mm / .022" x .028"

HSL 2723-13, Wire seizing and ligature pliers

Category: Pliers

Order number: HSL_2723-13



Weingart, sterilizable with TC-inserts, 135 mm / 5-3/8 inch Holds round, square and rectangular wires up to 0,56 mm x 0,70 mm / .022" x .028"

HSL 2720-13, Wire seizing and ligature pliers

Category: Pliers

Order number: HSL_2720-13



Weingart, sterilizable 130 mm (5 1/8 inch) Holds round, square and rectangular wires up to 0,56 mm x 0,70 mm / .022" x .028"

HSL 266-15, Separating pliers

Category: Pliers

Order number: HSL_266-15



150 mm

HSL 2226-14, Bracket removing pliers

Category: Pliers

Order number: HSL_2226-14




With TC-inserts 120 mm / 4-3/4 inch sterilizable

HSL 2225-99, Rubber replacement tip

Category: Pliers

Order number: HSL_2225-99



Rubber replacement tip for pliers HSL 2221-14. sterilizable, 1 Piece.

HSL 2225-14, Bracket removing pliers

Category: Pliers

Order number: HSL_2225-14



130 mm / 5-1/8 inch sterilizable

HSL 2221-14, Band removing pliers

Category: Pliers

Order number: HSL_2221-14



Oliver, sterilizable with TC-inserts 135 mm / 5-3/8 inch

HSL 222-14, Band removing pliers

Category: Pliers

Order number: HSL_222-14



OLIVER, no. 347, 145 mm 5-3/4

HSJ 155-04, Cheek retractor mirror

Category: Cheek / Lip retractors

Order number: HSJ_155-04



Fig. 4, 22 mm, plane (12 pcs,)

HSG 015-12 Probes

Category: Probes

Order number: HSG_015-12



Measuring probes, single ended, with 3 mm graduation, for the measurement of the orthodontic indication groups

HSC 508-15, Bracket placing tweezers

Category: Setting instruments

Order number: HSC_508-15



Hammacher, 130 mm (5 1/8 inch), spring action

HSC 507-15, Bracket placing tweezers

Category: Setting instruments

Order number: HSC_507-15



Bachmann, 140 mm (5 1/2 inch)

HSC 505-12, Bracket placing tweezers

Category: Setting instruments

Order number: HSC_505-12



145 mm (5 3/4 inch), spring action

HWC 500-15, Bracket placing tweezers

Category: Tweezers

Order number: HWC_500-15



Schrems, 150 mm, for placing molar brackets, from WIRONIT steel

HSC 509-12, Bracket placing tweezers

Category: Tweezers

Order number: HSC_509-12



125 mm (4 7/8 inch), indirect bonding tweezers

HSC 506-12, Bracket placing tweezers

Category: Tweezers

Order number: HSC_506-12



125 mm (4 7/8 inch)

HSC 504-13, Bracket placing tweezers

Category: Tweezers

Order number: HSC_504-13



130 mm (5 1/8 inch), double angled, crossed

HSC 503-00, Bracket placing tweezers

Category: Tweezers

Order number: HSC_503-00



130 mm (5 1/8 inch), crossed, self retaining / spring action

HSC 497-13, Bracket placing tweezers

Category: Tweezers

Order number: HSC_497-13



130 mm (5 1/8 inch), double angled, crossed, with cleaning scaler

HSJ 107-00, Mouth mirror handle

Category: Mouth mirrors and handles

Order number: HSJ_107-00



Mouth mirror measuring handle, hollow

HSL 252-20, Measuring ruler

Category: Measuring instruments

Order number: HSL_252-20



20 cm (8 inch)

HSL 252-15, Measuring ruler

Category: Measuring instruments

Order number: HSL_252-15



15 cm (6 inch)

HSL 248-70, Measuring ruler

Category: Measuring instruments

Order number: HSL_248-70



110 mm / 4-5/16 inch, stainless, sterilizable Measuring range: 0-60 mm

HSL 248-14, Space divider

Category: Measuring instruments

Order number: HSL_248-14



Korkhaus

HSL 247-52, Sliding caliper

Category: Measuring instruments

Order number: HSL_247-52



Zürcher Modell, 125 mm (4 3/4 inch)

HSL 247-51, Sliding caliper

Category: Measuring instruments

Order number: HSL_247-51



Münchner Modell, 125 mm (4 3/4 inch)

HSL 247-50, Sliding caliper

Category: Measuring instruments

Order number: HSL_247-50



Beerendonk, 125 mm (4 3/4 inch)

HWL 165-00, Ligature tucker

Category: Ligature instruments

Order number: HWL_165-00



Ligature tucker .024" / .016", from WIRONIT steel

HWL 164-00, Ligature tucker

Category: Ligature instruments

Order number: HWL_164-00



Ligature tucker .022" / .016", from WIRONIT steel

HSL 2743-14, Instrument for ligatures

Category: Ligature instruments

Order number: HSL_2743-14



Tying and elastomeric pliers, Mathieu, 140 mm, with spring, with very fine tips

HSL 2741-14, Instrument for ligatures

Category: Ligature instruments

Order number: HSL_2741-14



Tying and elastomeric pliers, Mathieu, 140 mm, with spring

HSL 264-06, Ligature tucker

Category: Ligature instruments

Order number: HSL_264-06



Ligature tucker

HSL 264-05, Ligature instrument

Category: Ligature instruments

Order number: HSL_264-05



for placing dentalastics

HSL 264-04, Ligature instrument

Category: Ligature instruments

Order number: HSL_264-04



Scaler

HSL 264-03, Ligature tucker

Category: Ligature instruments

Order number: HSL_264-03



with scaler

HSL 264-02, Ligature tucker

Category: Ligature instruments

Order number: HSL_264-02



HSL 264-01, Ligature tucker

Category: Ligature instruments

Order number: HSL_264-01



with band pusher

HSM 102-10 Sterilization Cassettes

Category: Cassettes / Boxes

Order number: HWM__102-10




Instrument holder with cover, stainless steel and silicon ,for 10 instruments

HSK 564-04, Punch

Category: Instruments for the invisible splint therapy

Order number: HSK_564-04



SLOT MACHINE acc. to Dr. Schwarze

HSK 561-01 Punch "The Gun 1"

Category: Instruments for the invisible splint therapy

Order number: HSK_561-01



Punch 200 mm length of the shank + 4 mm wide

HSB 630-10, Scissors for foils and splints

Category: Instruments for the invisible splint therapy

Order number: HSB_630-10



105 mm, curved, with micro saw edge, stainless

HWL 004-00, Inserting combination instrument

Category: Instruments for the lingual orthodontia

Order number: HWL_004-00



Inserting combination instrument

HWL 003-00, Inserting combination instrument

Category: Instruments for the lingual orthodontia

Order number: HWL_003-00




Inserting combination instrument

HSL 804-00, Lingual servant

Category: Instruments for the lingual orthodontia

Order number: HSL_804-00



Lingual servant acc. to Dr. Wilhelmy

HSL 803-00, Activating lever

Category: Instruments for the lingual orthodontia

Order number: HSL_803-00



Activating lever .016" x .022" acc. to Dr. Wilhelmy

HSL 802-00, Activating lever

Category: Instruments for the lingual orthodontia

Order number: HSL_802-00



Activating lever .014" acc. to Dr. Wilhelmy

HSL 801-00, Arch tucker

Category: Instruments for the lingual orthodontia

Order number: HSL_801-00



Arch tucker .025" acc. to Dr. Wilhelmy

HSL 4191-12A, Aderer

Category: Instruments for the lingual orthodontia

Order number: HSL_4191-12A



Wire and clasp bending pliers, 120 mm (4 3/4 inch), 90° angled, sterilizable For hard wire until 0,7 mm / .028" max. diameter

HSL 4191-12, Wire and clasp bending pliers

Category: Instruments for the lingual orthodontia

Order number: HSL_4191-12



Wire and clasp bending pliers, 120 mm / 4-3/4 inch, sterilizable, 90° angledSet from HSL 4191-12A and HSL 4191-12BFor hard wire until 0,7 mm / .028" max. diameter

HSL 4111-16B, Wire and clasp bending pliers

Category: Instruments for the lingual orthodontia

Order number: HSL_4111-16B



ADERER, 165 mm (6 1/2 inch), 90° angled, extra fine, sterilizable For hard wire until 0,4 mm / .016" diameter

HSL 4111-16A, Wire and clasp bending pliers

Category: Instruments for the lingual orthodontia

Order number: HSL_4111-16A



ADERER, 165 mm (6 1/2 inch), 90° angled, extra fine

HSL 4111-16, Aderer Set

Category: Instruments for the lingual orthodontia

Order number: HSL_4111-16



Wire and clasp bending pliers, extra fine , 90° angledSet from HSL 4111-16A and HSL 4111-16BFor hard wire until 0,4 mm / .016" diameter

HSL 3316-14, Distal end cutter

Category: Instruments for the lingual orthodontia

Order number: HSL_3316-14



With TC-inserts, 150 mm (5 7/8 inch) and safety hold. Cuts the wire directly behind the tubes or brackets - the cut wire is being held by the safety hold. -For braided wire from min. 0,36 mm / .014" Durchmesser until max. 0,4 x 0,4 mm / .016 x .016"-For round

HSL 3315-14, Distal end cutter

Category: Instruments for the lingual orthodontia

Order number: HSL_3315-14



Distal end cutter with safety hold, sterilizable With TC-inserts 140 mm / 5- 1/2 inch For hard wire from min. 0,46 mm / .018" until max 0,5 mm / .020" diameter For steel wire from min. 0,4 x 0,4 mm / .016 x .016" until max. 0,43 x 0,64 mm / .017 x .025"

HSL 2827-16, Ligature cutter

Category: Instruments for the lingual orthodontia

Order number: HSL_2827-16



These pliers are especially appropriate for the lingual technique. with TC-inserts, 165 mm (6 1/2 inch), with outside cutting-edge for soft wire from min. 0,25 mm / .010" until max. 0,36 mm / .014" diameter

HSL 2744-14, Instrument for ligatures

Category: Instruments for the lingual orthodontia

Order number: HSL_2744-14



Tying and elastomeric pliers, Mathieu, 140 mm, with TC-inserts, with spring, with very fine tips

HSJ 007-00, Orthodontic probe

Category: Instruments for the lingual orthodontia

Order number: HSJ_007-00




155 mm (6 1/8 inch), extra short

HSE 953-12, Instrument for ligatures

Category: Instruments for the lingual orthodontia

Order number: HSE_953-12



Mosquito clamp, modified, Wilhelmy/Hammacher, 120 mm, curved

HSE 950-12, Instrument for ligatures

Category: Instruments for the lingual orthodontia

Order number: HSE_950-12



Mosquito clamp, modified, Wilhelmy/Hammacher, 120 mm, straight

HSC 496-12 Placing tweezers

Category: Instruments for the lingual orthodontia

Order number: HSC_496-12




Placing tweezers, 125 mm, for an easy placing of the brackets- especially for bonding tubes

HSL 158-00, Band / bracket positioning gauge

Category: Band setting instruments

Order number: HSL_158-00



Band / bracket positioning gauge

HWL 389-18, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_389-18



UK-18P / XS - lower jaws / children size, perforated

HWL 388-17, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_388-17



UK-17P / S - lower jaws / children size, perforated

HWL 387-16, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_387-16



UK-16P / M - lower jaws / children size, perforated

HWL 386-15, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_386-15



UK-15P / L - lower jaws / children size, perforated

HWL 384-18, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_384-18



OK-18P / XS - upper jaws / children size, perforated

HWL 383-17, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_383-17



OK-17P / S - upper jaws / children size, perforated

HWL 382-16, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_382-16



OK-16P / M - upper jaws / children size, perforated

HWL 381-15, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_381-15



OK-15P / L - upper jaws / children size, perforated

HWL 379-18, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_379-18



UK-18 / XS - lower jaws / children size, regular

HWL 378-17, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_378-17



UK-17 / S - lower jaws / children size, regular

HWL 377-16, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_377-16



UK-16 / M - lower jaws / children size, regular

HWL 376-15, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_376-15



UK-15 / L - lower jaws / children size, regular

HWL 374-18, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_374-18



OK-18 / XS - upper jaws / children size, regular

HWL 373-17, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_373-17



OK-17 / S - upper jaws / children size, regular

HWL 372-16, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_372-16



OK-16 / M - upper jaws / children size, regular

HWL 371-15, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_371-15



OK-15 / L - upper jaws / children size, regular

HWL 334-99, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_334-99



99P, lateral, with swivel handle, partially toothed, perforated

HWL 333-32, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_333-32



32P UK, lower jaws, anterior, partially toothed, perforated

HWL 332-33, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_332-33



33P OK, upper jaws, anterior, partially toothed, perforated

HWL 331-31, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_331-31



31P right, partially toothed, perforated

HWL 330-30, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_330-30



30P left, partially toothed, perforated

HWL 325-99, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_325-99



99, lateral, with swivel handle, partially toothed, regular

HWL 324-32, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_324-32



32 UK, lower jaws, anterior, partially toothed, regular

HWL 323-33, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_323-33



33 OK, upper jaws, anterior, partially toothed, regular

HWL 322-31, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_322-31



31 right, partially toothed, regular

HWL 321-30, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_321-30



30 left, partially toothed, regular

HWL 319-04, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_319-04



UK-4P / S - lower jaws / perforated

HWL 318-03, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_318-03



UK-3P / M - lower jaws / perforated

HWL 317-02, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_317-02



UK-2P / L - lower jaws / perforated

HWL 316-01, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_316-01



UK-1P / XL - lower jaws / perforated

HWL 315-05, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_315-05



UK-5 / XS - lower jaws / regular

HWL 314-04, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_314-04



UK-4 / S - lower jaws / regular

HWL 313-03, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_313-03



UK-3 / M - lower jaws / regular

HWL 312-02, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_312-02



UK-2 / L - lower jaws / regular

HWL 311-01, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_311-01



UK-1 / XL - lower jaws / regular

HWL 309-04, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_309-04



OK-4P / S - upper jaws / perforated

HWL 308-03, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_308-03



OK-3P / M - upper jaws / perforated

HWL 307-02, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_307-02



OK-2P / L - upper jaws / perforated

HWL 306-01, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_306-01



OK-1P / XL - upper jaws / perforated

HWL 305-05, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_305-05



OK-5 / XS - upper jaws / regular

HWL 304-04, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_304-04



OK-4 / S - upper jaws / regular

HWL 303-03, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_303-03



OK-3 / M - upper jaws / regular

HWL 302-02, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_302-02



OK-2 / L - upper jaws / regular

HWL 301-01, Wiro-lock impression trays

Category: Diagnosis

Order number: HWL_301-01



OK-1 / XL - upper jaws / regular

HSL 248-50, Interstice width caliper

Category: Diagnosis

Order number: HSL_248-50



Measuring range: 4 - 80 mm

По вопросам продаж и поддержки обращайтесь:

Алматы (727)345-47-04
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(727)345-47-04

Беларусь +(375)257-127-884

Узбекистан +998(71)205-18-59

Киргизия +996(312)96-26-47

эл.почта: hqs@nt-rt.ru || сайт: <https://hammacher.nt-rt.ru/>